

# MRI Scan Patient Checklist

## remedē® System

Patient name: \_\_\_\_\_ Patient date of birth: \_\_\_\_\_

Today's date (pre-MRI check): \_\_\_\_\_

Scheduled MRI date: \_\_\_\_\_

### Part A – Checklist for the pre-MRI check

	YES	NO
1. The patient was told to bring their patient ID card to the MRI scan appointment.	<input type="checkbox"/>	<input type="checkbox"/>
2. Confirmed the patient's implanted <b>remedē</b> ® System components are MR Conditional (see table 1 and table 2 on the back of this page).	<input type="checkbox"/>	<input type="checkbox"/>
3. Confirmed the patient doesn't have any abandoned respistim leads.	<input type="checkbox"/>	<input type="checkbox"/>
4. Confirmed whether or not patient has other implantable device(s). For non-ZOLL Respicardia implantable devices and components, a healthcare provider will need to ensure MRI compatibility.	<input type="checkbox"/>	<input type="checkbox"/>
5. The <b>remedē</b> System was interrogated and the patient's therapy settings were stored (this information may be used for restoring the patient's therapy following the MRI scan).	<input type="checkbox"/>	<input type="checkbox"/>
6. A lead integrity interrogation was performed using the Physician Programmer.	<input type="checkbox"/>	<input type="checkbox"/>
7. All electrode impedances are between the <b>remedē</b> System Min/Max impedance limits.	<input type="checkbox"/>	<input type="checkbox"/>

If all the answers above are **YES**, complete Part B and scan the form into your Electronic Medical/Health Record system and if necessary, fax it to the imaging center performing the patient's MRI scan.

If any answer is **NO**, do not proceed with the MRI scan for this patient.

### Part B - Healthcare provider's acknowledgment of MRI scan readiness for the remedē System and components

I \_\_\_\_\_ confirm that the **remedē** System has been checked prior to the MRI scan.

- ☐ The patient doesn't have any abandoned respistim leads.
- ☐ At the time of the check, all electrode impedances for the **remedē** System were within Min/Max limits.

**Note:** The **remedē** System may enter into Safe Mode and display an alert in the presence of a very large magnetic field such as an MRI machine. Please contact your ZOLL Respicardia representative.

**Table 1**

**remedē System components that are eligible for Whole-body MRI scans and Head Scan (1.5T and 3T) under specified conditions.**

COMPONENT	MODEL NUMBER(S)
remedē IPG	1001, 1100, 1600
respistim Left Stimulation Lead	L Lead Models – 2002, 2003, 2004 LQ Lead Models – 50XX, 51XX, 56XX LQS Lead Models – 40XX, 41XX, 46XX Where “XX” denotes Lead Length in cm (45, 55, 65 or 85 cm)
respistim Right Stimulation Lead	R Lead Models - 3101–3106, 3201–3206, 3601–3606, 3611–3616, 3651–3656

**Table 2**

**Off-the-shelf components that are eligible for Whole-body and Head MRI scans (1.5T and 3T) when used with the remedē System under specified conditions.**

COMPONENT DESCRIPTION	USE CONSIDERATIONS
Bipolar Sensing Lead with IS-1 terminal and co-radial coil	The <b>remedē</b> Model 1001/1100 IPGs were evaluated to be MR Conditionally safe using an off-the-shelf, IS-1 bipolar lead with co-radial coil connected to the IPG sensing port. However, prior to an MRI examination, determine whether the patient has a sensing lead inserted into the <b>remedē</b> Model 1001/1100 IPG sensing port. If so, the MR labeling of the sensing lead must be determined and the most restrictive MRI exposure requirements must be used of the medical device implants. Contact the appropriate device manufacturers if you have questions. If you are unclear what implants may be present, perform an x-ray to determine implant type and location. Do not conduct an MRI examination if any conditions or implants that would prohibit or contraindicate an MRI are present.
IS-1 Port Plug (3.5cm) with metallic contacts	The <b>remedē</b> Model 1001 / 1100 IPGs were evaluated to be MR Conditionally safe using an off-the-shelf, IS-1 port plug (3.5 cm in length) with metallic contacts inserted into the IPG sensing port. However, prior to an MRI examination, determine whether the patient has an IS-1 port plug inserted into the <b>remedē</b> Model 1001 / 1100 IPG sensing port. If so, ensure the IS-1 port plug is no greater than 3.5cm in length or has MR conditional approval. Contact the appropriate device manufacturers if you have questions. If you are unclear what implants may be present, perform an x-ray to determine implant type and location. Do not conduct an MRI examination if any conditions or implants that would prohibit or contraindicate an MRI are present.

#### Additional information regarding MRI guidelines for the remedē System



Visit the **remedē** website to access the full **remedē MRI Guidelines Manual** or use the QR code

For any questions, call ZOLL Respicardia at 952-540-4470.

## ZOLL Respicardia, Inc.

12400 Whitewater Dr., Suite 150 | Minnetonka, MN 55343  
952-540-4470 | [info@remede.zoll.com](mailto:info@remede.zoll.com) | [remede.zoll.com](http://remede.zoll.com)

Copyright © 2025 ZOLL Medical Corporation. All rights reserved. Respicardia and remedē are registered trademarks of ZOLL Respicardia, Inc. in the United States and/or other countries. ZOLL is a registered trademark of ZOLL Medical Corporation in the United States and/or other countries. MedEd3147, Rev B

# ZOLL®